



Gentueri Oral Collection Kit

The Gentueri Oral Collection Kit is designed for the simple collection, preservation and processing of samples for genetic testing. Your sample will be protected and preserved at room temperature without any messy liquids. This is the easiest non-invasive collection kit to use and process.

The kits are designed to simplify the workflow by allowing you to eject the swab head of the SimpleSwab™ into the SwabSaver® directly once the sample has been collected. At this point, the SwabSaver® is like a 2-milliliter, micro-centrifuge tube with a spin basket and DNA extraction can be performed without transferring the sample. This will save time and reduce the potential for transfer errors in the laboratory.



Kit Contents	Description
SwabSaver®	Protects and preserves swab at room temperature
SimpleSwab™	Sterile, DNA-free swab with ejectable head
Instructions	Easy to follow directions to properly use the kit
Return Pouch	Small return pouch reduces shipping costs

Our unique kit design preserves samples at room temperature and allows plenty of sample to be collected with one swab. Shipping is less costly in a small return pouch without liquid preservative. Sample tracking is easy with matching barcodes printed on the SwabSaver, kit box, return pouch and on the instructions kept by the user to reference for obtaining results.

The Gentueri Oral Collection Kit is customizable to individual needs. We have experience customizing many types of kits, from sexual assault sample collection kits to kits from collecting samples for genetic testing. Let us know what customized features you require to simplify your workflow.

Item	Quantity	Catalog Number
Gentueri Oral Collection Kit	25/CS	5000-025
Gentueri Oral Collection Kit	100/CS	5000-100
Gentueri Oral Collection Kit	XXX/CS (Custom)	5000-XXX

Catalog and order information: www.Gentueri.com
 Phone: 608-838-1640 Email: info@gentueri.com