



twitter com/gentueriing

GenSwab[™] in partnership with → AHLSTROM MUNKSJÖ



- Easily collect an oral sample in 1 minute
- Biobank at ambient temperature for over 20 years
- Sterile swab with color change indicator on the card
- Dry preservation by desiccation on a solid support
- Barcoding for a chain of custody
- Sterilized and manufactured in the USA

Complementary Products

GenSwab[™] Variant

GenSwab

GenSwab 2

- GenDry[™] GS pouch (Cat# G-7501-000) for transport from the collection site to the laboratory
- GenPouch[™] Large, Heat Seal Closure (Cat# G-7201-000) for long-term storage of processed specimen cards

Card Size

3" x 2"

2" x 2"

- GenPouch[™] Large, Tape Closure (Cat# G-7203-000) for long-term storage of processed specimen cards
- GenPouch[™] Small, Tape Closure (Cat# G-7204-000) for long-term storage of processed specimen cards
- GenKeep $^{\scriptscriptstyle\rm TM}$ (Cat# G-9107-000) for cost-effective long-term storage of specimen cards

Catalog #

G-7300-100

G-7400-100

• Customizable mailer and return envelope (Cat# N/A)



+1-608-949-8568 info@gentueri.com www.gentueri.com

